

## Protease Enzymes and Jello


### Materials

- 1 packet of Jello powder
- Deionized water
- Weigh Boats
- 250-mL Erlenmyer Flask
- 10-mL Graduated Cylinder
- Pasteur Pipettes and bulbs
- Spectrophotometer test tubes
- Spectrophotomer
- Heating Block
- Papain enzyme solution

### Procedure

1. Jello Preparation
  - a. Stir one packet of Jello into 350 mL of hot deionized water
  - b. When Jello is completely dissolved, allow to cool for 10 minutes
  - c. Pour Jello into weigh boats approximately 1 cm deep
  - d. Allow to thicken and dry overnight
2. Jello Absorption Spectrum
  - a. Use a moistened test tube to cut out one piece of Jello (as pictured below)
  - b. Place the piece into a 250-mL Erlenmeyer flask
  - c. Add 3mL deionized water
  - d. Swirl the contents for 3 minutes
  - e. Zero the spectrophotometer at wavelength=460 nm as described in the standard operating procedure
  - f. Pipette the liquid into a clean test tube being careful not draw up Jello particles (having particles in solution will cause error in absorbance measurements)
  - g. Insert the tube in the spectrophotometer. Measure and record the absorbance in table 1 in the data sheets
  - h. Repeat steps e.-g. for wavelengths 470 nm, 480 nm, 490 nm, 500 nm, 510 nm, 520 nm, and 530 nm using the same test tube of solution
  - i. Make a graph of Absorbance vs. wavelength and determine the wavelength of maximal absorbance ( $\lambda_{\text{max}}$ )
3. Spectrophotometer Measurements
  - a. Use a moistened test tube to cut out one piece of Jello
  - b. Place into 250-mL Erlenmeyer flask
  - c. Add 3 mL deionized water
  - d. Swirl contents for 1 minute timed and record observations
  - e. Pipette the liquid into a test tube, being careful not to draw up Jello particles (as shown below)
  - f. Zero the spectrophotometer for the  $\lambda_{\text{max}}$  as described in the standard operating procedure
  - g. Insert the test tube into the spectrophotometer. Measure and record its absorbance
  - h. Remove the test tube and pour the solution back into the Erlenmeyer flask
  - i. Add 3 mL of the Papain solution to the Erlenmeyer

- j. Repeat steps d.-h. 10 times, recording all measurements and observations in table 2 in the data section
  - k. Make a graph of Absorbance vs. time swirled and describe the trend
4. Measurements with Heat
- a. Use a moistened test tube to cut out another piece of Jello
  - b. Pour 3 mL of deionized water into the test tube
  - c. Place in the heating block and wait for 1 minute
  - d. Pipette the solution into another test tube, being careful once again to not take up any Jello particles
  - e. Be sure the spectrophotometer is zeroed
  - f. Measure the absorbance of the solution
  - g. Return the solution to the test tube with Jello
  - h. Repeat steps c.-g. 10 times, recording measurements and observations in table 3
  - i. Make a graph of Absorbance vs. time heated and describe the trend



## Data Sheet

Table 1 Absorption Spectrum

Wavelength (nm)	Absorbance
460	
470	
480	
490	
500	
510	
520	
530	

$\lambda_{\max}$ : \_\_\_\_\_

Table 2 Spectrophotometer Measurements:

<b>Measurement #</b>	<b>Time swirled (min:sec)</b>	<b>Absorbance</b>	<b>Appearance</b>
0 (no Papain)			
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

Data trend:

Table 3 Measurements with Heating

<b>Measurement #</b>	<b>Time heated (min:sec)</b>	<b>Absorbance</b>	<b>Appearance</b>
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			


Data Trend:

### Jello and Papain Supplement

Here is some example data from when I ran the experiment:

Part 1)

Wavelength (nm)	Absorbance
460	0.114
470	0.102
480	0.151
490	0.259
500	0.272
510	0.229
520	0.164
530	0.145




$\lambda_{\max}=500 \text{ nm}$

Part 2)

Measurement #	Time swirled (min:sec)	Absorbance	Appearance
0 (no Papain)	1:00	0.014	
1	1:00	0.124	


2	2:00	0.153	
3	3:00	0.154	
4	4:00	0.178	
5	5:00	0.22	
6	6:05	0.257	
7	7:05	0.26	
8	8:05	0.391	
9	9:05	0.402	
10	10:05	0.468	



Part 3)

Measurement #	Time heated (min:sec)	Absorbance	Appearance
1	1:05	0.049	
2	2:06	0.125	
3	3:06	0.166	
4	4:12	0.349	
5	5:12	0.398	
6	6:15	0.535	
7	7:15	0.599	

8	8:15	0.596	
9	9:15	0.696	
10	10:15	0.655	



**Note:** Another variation of part 3 could be to use an acid or base to change the pH and cause protein or enzyme denaturation. We did not investigate this yet because of time constraints.

#### **VWR Product List**

##### *Package 6 (Papain):*

10 mL graduated cylinders	65000-000
13x100 mm test tubes	47729-572
PTFE stir rods	89026-280
Papain carica papaya, 50 g	80057-240

This experiment was adapted from:

Hagar, William G. and Bullerwell, Lornie D. Supermarket Proteases. *The Science Teacher*. October, 2003 pp. 26-30.